



ARBEITSAUFTRAG PATIENT ANLEGEN

1. Lege dich im Archiv als neuen Patienten an (Vorname, Nachname)
2. Lege im Archiv den Kunden an (Klasse, Lehrer)
3. Markiere im Zahnschema die Zähne welche gescannt werden müssen
4. Gib im Archiv den ausführenden Techniker ein (Lehrer)
5. Speichere das Datenblatt
6. Öffne die Scan - Software

ORDNE DEN FELDERN DIE ENTSPRECHENDEN ZIFFERN ZU!!

The screenshot shows the Zirkonzahn Archiv software interface. It features a top navigation bar with icons for 'Neu', 'Laden', 'Upload', 'Duplizieren', 'Scan', and 'Modeller'. Below this is a sidebar with sections for 'Kunde', 'Patient', and 'Techniker', each containing input fields for name, number, and other details. The main area displays a dental arch diagram with numbered teeth (1-32) and a 'Bitte zuerst einen Kunden auswählen' (Please first select a customer) message. Orange boxes and arrows are placed over the interface to indicate where to enter data: one box points to the 'Kunde' name field, another to the 'Patient' name field, and a third to the 'Techniker' name field. A fourth box points to the 'Scan' button in the top navigation bar. A fifth box points to the 'Zirkonzahn' logo in the bottom right corner of the interface.